



Graver Technologies

Filtration | Separation | Purification

QMA™ Polypropylene Filter Series

“Absolute” Rated High Performance Pleated Polypropylene Filter Cartridge

This filter is constructed with a high surface area melt blown polypropylene media for low initial pressure drop, high dirt holding capacity, and high efficiency performance.

Filter Features–Benefits

- Micron ratings from 0.2 to 20µm – Broad application range
- Meets current USP Class VI biological test for plastics
- “Absolute” Efficiency – Rated at 99.98% (Beta 5000)
- High surface area – High flow rate, and long service life – Minimize maintenance cost
- Fixed Pore construction – Resists dirt unloading at maximum differential pressure
- Polypropylene Construction – Inert to many process fluids
- Various Gasket/O-ring materials – Compatible with many fluids
- Heavy duty molded cage – High structural strength
- Highly consistent melt blown media for consistent performance

Filter Specifications

Construction material:	Polypropylene
Gaskets/O-Rings:	Buna-N, EPR, Silicone, Viton, Teflon Encapsulated Viton (O-Rings only)
Micron ratings:	0.2, 0.45, 1.0, 2.5, 5.0, 10, 20µm
Dimensions and Operating Parameters	
Nominal lengths:	9.75" 10", 20", 30", 40" (24.7, 25.4, 50.8, 76.2, 101.6 cm)
Outside diameter:	2.7" (6.86 cm)
Inside diameter:	1.0" (2.54 cm)
Surface area:	up to 7.0 ft ²
Maximum operating temperature:	176°F (80°C)
Maximum differential pressure:	75 psid @ 70°F (5.2 bar @ 21°C) 40 psid @ 176°F (2.8 bar @ 80°C)

FDA & USP compliance

All polypropylene material used in manufacturing complies with the regulations of the Food and Drug Administration (FDA) title 21 of the Code of Federal Regulations Sections 174.5, 177.1520, and 177.1630, as applicable for food and beverage contact.



Filter Removal Efficiency

Beta Ratio Efficiency	Beta 5000 99.98%	Beta 100 99%	Beta 50 98%
0.2 micron	0.20	0.10	0.05
0.45 micron	0.45	0.30	0.20
1.0 micron	1.0	0.60	0.30
2.5 microns	2.5	2.0	1.5
5.0 microns	5.0	4.0	3.0
10.0 microns	10.0	8.0	7.0
20.0 microns	20.0	17.0	15.0

$$\text{Beta Ratio} = \frac{\text{Upstream particle counts}}{\text{Downstream particle counts}}$$

The micron ratings shown at various efficiency and beta ratio value levels were determined through laboratory testing, and can be used as a guide for selecting cartridges and estimating their performance. Under actual field conditions, results may vary somewhat from the values shown due to the variability of filtration parameters.

Testing was conducted using the single-pass test method, water at 2.5 gpm/10" cartridge. Contaminants included latex beads, coarse and fine test dust. Removal efficiencies were determined using dual laser source particle counters.

Applications

- Food & beverage
- Bottled water
- Process water
- Aqueous solutions
- Pharmaceutical
- RO Prefilters
- Chemicals
- Cosmetics
- Ink

QMA Nomenclature Information


QMA	1	-20	P3	V	-R	-I
Filter Type QMA Series Filters		Nominal Length (inches) -9.75 -10 -20 -30 -40			Pre-Rinse Option -R Factory Pre-Rinse	Insert -I End cap insert for steaming
Retention Rating (microns) 0.2 0.45 1 2.5 5 10 20						
End Configuration P Double Open End P2 226/Flat Single Open End P3 222/Flat Single Open End P7 226/Fin Single Open End P8 222/Fin Single Open End PX Extended Core AM Single open end, internal O-Ring NPC Double open end, internal O-Ring			Gasket or O-Ring S Silicone B Buna-N E EPDM V Viton T Teflon encap. Viton (O-Rings only) Teflon (Gasket only)			

Example: QMA 1-20 P3V-R-I

Performance Specifications

Sanitization

Hot water at 176°F (80°C) at 5 psid (0.35 bar) for 30 min.
 In-line steam at 257°F (125°C) at 1 psid (0.07 bar) for 30 min.
 Autoclavable at 257°F (125°C) for 30 min.



For more information

Graver Technologies Customer Service: **1-888-353-0303**

Technical Support: **1-800-510-0932**

E-mail us at info@gravertech.com

Graver Technologies Europe (UK): **+44-1424-777791**

All information and recommendations appearing in this bulletin concerning the use of products described herein are based on tests believed to be reliable. However, it is the user's responsibility to determine the suitability for his own use of such products. Since the actual use by others is beyond our control, no guarantee, expressed or implied, is made by Graver Technologies as to the effects of such use or the results to be obtained. Graver Technologies assumes no liability arising out of the use by others of such products. Nor is the information herein to be construed as absolutely complete, since additional information may be necessary or desirable when particular or exceptional conditions or circumstances exist or because of applicable laws or government regulations.

QMA is a trademark of Graver Technologies, LLC.

DISTRIBUTED BY:



200 Lake Drive
 Glasgow,
 DE 19702 U.S.A.

302-731-1700
 800-249-1990
 Fax: 302-369-0938

e-mail: info@gravertech.com
 web site: www.gravertech.com

